

NEW

wiscan®

HERMES 24-7

Fully Autonomous Imaging System

The freedom to focus on the fun stuff



IDEA Bio-Medical Ltd.
Seeing Better


WISCAN®

HERMES 24-7 Fully Autonomous Imaging System

Minimal Involvement


With push-button programming the system autonomously:

- Incubates plates in a built-in plate stacker which accommodates up to 20 plates under incubation conditions (37°C, 5% CO₂)
- Loads & unloads plates for imaging
All external sources supported –Incubators, stackers, liquid handlers
- Acquires high quality images on the WiScan® Hermes automated high content screening system
- Analyzes acquired images using WiSoft® Athena automated image analysis software
- Produces presentation-ready results using multiple visualization tools (scatter plots, histograms, heat maps, time graphs, dose curves, etc.)



Compact

- Minimal footprint design
- Flexible access from any direction
- Reliable re-positioning




Hermes HCS
Imaging System

Dedicated Cart

Manipulator Arm

Precision
Docking Station



www.idea-bio.com

Worldwide headquarter

IDEA Bio-Medical Ltd.
2 Prof. Bergman St.
Rehovot, 76705, Israel
Phone: +972 732732 400
Fax: +972 89469 556
info@idea-bio.com

Asia - Pacific

Chen-Tech Challenge (HK) Limite
22/F Hing Yip Commercial CTR
280 Des Voeux Rd. Central, HK
Phone: +86 20 37651309
Maor-T@chen-tech.hk

USA

IDEA Bio-Medical Ltd.
5632 Van Nuys Blvd.
Suite # 411, Van Nuys, CA 91401
phone: 818-823-9763
info@idea-bio.com

NEW